VASCO

TECHNICAL DATA SHEET

All Clear Hayne Set

1. Description

The ALL CLEAR HAYNE Kit includes:
HAYNE Lens Cleaner (30 ml),
HAYNE MIST Anti-Fog Lens Cleaner (15 ml),
Microfibre cloth (180 × 150 mm)

2. Technical Data

1. Medical Device || HAYNE Lens Cleaner (30 ml) professional cleaning fluid for prescription glasses lenses. Also suitable for other optical surfaces. Protects against dirt and gently cares for the cleaned coatings. Alcohol-free.

Instructions for use: Spray HAYNE Lens Cleaner on the surface to be cleaned and wipe with a soft microfiber cloth. Do not use paper towels or tissues – this may damage the coating or scratch the lens. Avoid contact with eyes. If contact occurs, immediately rinse with plenty of water. Keep out of the reach of children.

2. HAYNE MIST Anti-Fog Lens Cleaner (15 ml) – a specialized spray that professionally protects all optical surfaces from fogging and additionally protects against dirt.

3. Microfiber Cloth (180 x 150 mm) – thanks to its unique structure, it cleans perfectly, eliminating the problem of stubborn streaks.

4. ALL CLEAR HAYNE Box – a special box measuring (H x W x L) 3 x 15 x 14 cm, made of white, coated paper, with three dedicated compartments for the above-mentioned products.

3. Producer

Manufacturer:
Hayne Polska Sp. z o.o.,
ul. Jana Henryka Dąbrowskiego 286/288,
60-406 Poznań

Importer:
Vasco Electronics S.A.
Al. 29-Listopada,
31-401 Krakow, Poland

4. Notes

This technical data sheet has been carefully drawn up to the best of our knowledge. The information provided is based on practical experience, in-house testing and reflects our current level of knowledge.
It is intended for information only and does not constitute a guarantee in terms of product properties or its suitability for specific applications. We accept no liability for any mistakes, errors in standards, or printing errors. The intellectual property rights of this document belong to our company. Any use beyond the authorized scope requires prior consultation. Unauthorized copying or reproduction is strictly prohibited.
This is a medical device. It should be used in accordance with the operating instructions or label.